



Specifications

Bag	Fill volume (average)	Admixture volume (average)	Residual volume (average)	Empty primary bag alone (average weight)	Empty primary bag with overwrap (average weight)
	ml	ml	ml	g	g
50 ml	59	86	0,9	8,1	13,4
100 ml	111	77	1,0	8,5	14,1
250 ml	271	154	1,3	11,3	19,0
500 ml	530	241	2,5	15,3	25,7
1000 ml	1047	172	4,4	18,6	30,3

Baxter data on file available on request

Baxter Medical AB
Box 63
164 94 Kista

www.baxter.se

Baxter A/S
Tobaksvejen 2A
2860 Søborg

www.baxter.dk

Baxter Oy
PL 119
00181 Helsinki

www.baxter.fi

Baxter AS
Gjerdrumsvei 11
0484 Oslo

www.baxter.no

Baxter